

TO BE A WORLD-LEADING ANALYTICAL TESTING SOLUTIONS PROVIDER

- Spectroscopy
- Chromatography
- Mass Spectrometry



EDX2800B
EDX Spectrometer

Instrument Introduction:

With the widespread of EDX2800B over different fields, we design this type to meet the need of optimizing the product performance and improve the safety protection grade.


The reliability of the product is improved by using the high voltage source and copper tube of the new generation and the testing efficiency is improved by the adopting the high power



Testing 75 kinds of elements, 1ppm limit of detection, Repeatability 0.1%, Stability 0.1%



ROHS



Application field :

- RoHS testing

The performance is stable and reliable, achieving higher test efficiency

The automatic function of door sensor and high voltage lock gives you protection from all directions

Performance advantages:

- Down-side copperSource: meet the test requirements of samples of different kinds and shapes
- Collimator and filter: the Auto-switch between various collimators and filters to meet the application of different testing methods
- Movable platform: sophisticated manual movable platform is convenient for locating test point
- High-resolution detector: improve the analyzing accuracy
- High voltage source and copper tube of the new generation: the performance is stable and reliable, achieving higher test efficiency

Technical specifications:

- Measurable elements: S to U.
- Limit of detection (LOD) reaches 1ppm.
- Element content: ppm to 99.99%
- Arbitrary optional analysis and identification models
- Independent matrix effect correction models
- Multi-variable non-linear regression procedure

Standard configurations:

- Movable sample platform
- Signal-to-Noise Enhancer (SNE)
- Electric-cooling Si-PIN detector
- Signal detection electronic circuit
- High and low voltage power

- Ambient temperature: 15°C to 30°C
- Power supply: AC 220V±5V, AC purified stabilized voltage power supply.
- Energy resolution: 160±5eV
- Sample chamber size: 439mm×300mm×50mm
- Instrument size: 550mm×410mm×320mm
- Instrument weight: 45kg

Streamlined man-machine design, promising your operation security


Operation indicator makes you operate comfortably



Test cases:




RoHS Testing



Test Results Spectrum



Plating Thickness Testing



Test Results Spectrum

What are RoHS and WEEE Directives?

The European Union has adopted Directive 2002/95/EC on the restriction of certain hazardous substances (RoHS) and Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) with their publication in the Official Journal of the European Union on February 13, 2003. WEEE comes into effect on August 13, 2005 and RoHS requires the substitution of various heavy metals (lead, mercury, cadmium and hexavalent chromium) and brominated flame retardants (polybrominated biphenyls [PBB] or polybrominated diphenyl ethers [PBDE]) in new electrical and electronic equipment put on the market from July 1, 2006.



Testing standard of substances restricted by RoHS Directive

Hazardous substances	Standards (mg/kg)
Cd	100
Pb	1000
Hg	1000
Cr ⁶⁺	1000
PBBs	1000
PBDEs	1000

Restricted substances and their typical uses

Pb	
Solders	
Paints	Pigments and driers
Glass materials	Pb is allowed in fluorescent lamp
Ceramic materials	Pb is allowed in certain electronic ceramic materials
Iron, aluminum and copper materials	A certain amount of Pb is allowed
Plastics	PVC stabilizer and pigments
Batteries	Pb is allowed in acidic batteries for vehicles
Cd	
Plastics	Stabilizer and pigments
Solders	Seldom used
Ceramics	Seldom used
Connectors	Relays and switches
Batteries	Cd is allowed in Ni-Cd batteries
Semiconductors	Optical sensors and solar cell panels
Hg	
Batteries	Prohibited (see battery directive)
Connectors	Relays and sensitive switches
Fluorescent lamps	A certain amount of Hg is allowed
Cr ⁶⁺	
Passivation layers	Commonly used for naked metal surfaces to enhance adhesion of plating layers
Anti-corrosive plating layers	Painting and plating layers
Chrome plating layers	Plating layer of chromium metal is not under control
Plasticizer	Commonly used to plastics plating process but not final products
PBBs & PBDEs	
Plastics	Brominated flame retardants

The analytical method of filtration for XRF to test RoHS substances



PERIODIC TABLE OF ELEMENTS

Per	IA																	0	
1	1 H 1.008																	2 He 4.008	
2	3 Li 6.94 0.052	4 Be 9.012 0.110											5 B 10.81 0.185	6 C 12.01 0.282	7 N 14.01 0.392	8 O 15.99 0.523	9 F 18.99 0.677	10 Ne 20.17 0.861	
3	11 Na 22.99 1.041	12 Mg 24.31 1.254			III B	IV B	V B	VI B	VII B	VIII		I B	II B	13 Al 26.99 1.487	14 Si 28.09 1.740	15 P 30.97 2.015	16 S 32.06 2.307	17 Cl 35.45 2.622	18 Ar 39.94 2.957
4	19 K 39.1 3.312 3.589	20 Ca 40.08 3.690 4.012	21 Sc 44.96 4.088 4.459	22 Ti 47.90 4.508 4.931	23 V 50.94 4.949 5.427	24 Cr 51.99 5.411 5.947	25 Mn 54.94 5.895 6.492	26 Fe 55.84 6.400 7.059	27 Co 58.93 6.925 7.649	28 Ni 58.7 7.472 8.265	29 Cu 63.54 8.041 8.907	30 Zn 65.38 8.631 9.572	31 Ga 69.72 9.243 10.26	32 Ge 72.5 9.876 10.98	33 As 74.92 10.53 11.73	34 Se 78.9 11.21 12.50	35 Br 79.90 11.91 13.29	36 Kr 83.8 12.63 14.12	
5	37 Rb 85.47 13.38 14.97 1.694 1.752	38 Sr 87.82 14.14 15.85 1.806 1.872	39 Y 88.91 14.93 16.75 1.922 1.996	40 Zr 91.22 15.75 17.69 2.124 2.257	41 Nb 92.91 16.58 18.65 2.266 2.402	42 Mo 95.94 17.44 19.63 2.293 2.424	43 Tc #(99)	44 Ru 101.0 20.17 21.69 2.558 2.696	45 Rh 102.9 21.12 22.76 2.696 2.838	46 Pd 106.4 22.10 23.86 2.984 3.133	47 Ag 107.9 23.11 24.99 3.133 3.287	48 Cd 112.4 23.11 26.14 3.133 3.287	49 In 114.8 24.14 27.38 3.92 4.131	50 Sn 118.6 25.19 28.60 4.131 4.347	51 Sb 121.7 26.27 29.85 4.347 4.570	52 Te 127.6 27.38 31.13 4.570 4.800	53 I 126.9 28.51 32.44 4.800 5.036	54 Xe 131.3 29.67 33.78 5.036 5.366	
6	55 Cs 137.3 30.85 35.15 4.286 4.620 5.280 3.794	56 Ba 137.3 32.07 36.55 4.467 4.828 5.531 3.953	Ln		72 Hf 178.4 55.38 63.56 7.898 9.021 10.51 6.958	73 Ta 180.9 57.11 65.56 8.145 9.341 10.81 7.172	74 W 183.8 58.86 67.59 8.396 9.670 11.28 7.386	75 Re 186.2 60.66 69.66 8.651 10.01 11.68 7.602	76 Os 190.2 62.48 71.78 8.910 10.35 12.09 7.821	77 Ir 192.2 64.35 73.93 9.173 10.71 12.51 8.040	78 Pt 195.0 66.25 76.13 9.441 11.07 13.38 8.267	79 Au 197.0 68.19 78.37 9.711 11.44 14.28 8.493	80 Hg 200.5 70.16 80.66 9.987 11.82 15.24 8.720	81 Tl 204.3 72.18 82.99 10.27 12.61 14.76 8.952	82 Pb 207.2 74.23 85.36 10.55 13.02 15.74 9.183	83 Bi 208.9 76.32 87.77 10.84 13.02 15.74 9.419	84 Po #(209)	85 At #(210)	86 Rn #(222)
7	87 Fr #(223)	88 Ra #(226)	An																

- Alkali Metals
- Non-Metal
- Transitional element
- Halogen
- Lanthanoids
- Actinides
- Alkaline earth
- Rare gases
- Main group metal

Mark 1: #Radioactive Elements *Man Made Elements
 2: All the Numbers are ordered one by one in this way, Atomic Number
 Element Symbol, Atomic Weight, α, Kβ, Lα, Lβ, Lγ, Le

Ln	57 La	58 Ce	59 Pr	60 Nd	61 Pm	62 Sm	63 Eu	64 Gd	65 Tb	66 Dy	67 Ho	68 Er	69 Tm	70 Yb	71 Lu
	138.9	140.1	140.9	144.2	#(147)	150.4	152.0	157.2	158.9	162.5	164.9	167.2	168.9	173.0	175.0
	33.30	34.57	35.86	37.19	38.54	39.91	41.32	42.76	44.23	45.73	47.26	48.82	50.41	52.04	53.59
	37.99	39.45	40.95	42.48	44.05	45.65	47.28	48.95	50.65	52.38	54.16	55.96	57.81	59.69	61.61
	4.651	4.840	5.034	5.230	5.431	5.636	5.846	6.059	6.275	6.495	6.720	6.948	7.181	7.414	7.654
	5.043	5.262	5.489	5.722	5.956	6.206	6.456	6.714	6.979	7.249	7.528	7.810	8.103	8.401	8.708
5.789	6.052	6.322	6.602	6.891	7.180	7.478	7.778	8.104	8.418	8.748	9.089	9.424	9.779	10.14	
4.124	4.287	4.452	4.632	4.816	4.994	5.176	5.361	5.546	5.742	5.942	6.152	6.341	6.544	6.752	
An	89 Ac	90 Th	91 Pa	92 U	93 Np	94 Pu	95 Am	96 Cm	97 Bk	98 Cf	99 Es	100 Fm			
	#(227)	(232)	231.0	238.0	(237)	(244)	(243)	(247)	(247)	(251)	(252)	(257)			
	89.79	92.19	94.64	97.14	99.69	102.3	104.9	107.7	110.5	113.3	116.2	119.2			
	103.3	106.1	108.9	111.8	114.7	117.7	120.8	123.9	127.1	130.4	133.7	137.2			
	12.65	12.97	13.29	13.61	13.95	14.28	14.62	14.96	15.31	15.66	16.02	16.38			
	15.71	16.2	16.7	17.22	17.74	18.28	18.83	19.39	19.97	20.56	21.17	21.79			
18.41	18.98	19.55	20.16	20.77	21.40	22.04	22.69	23.37	24.06	24.76	25.47				