

SK-L330-Pro

Features:

- ♦ Max. load capacity is 7.5kg(L330) and 2.5kg(L180) with platform.
- ◆ Double LCDs for independent display of time and speed.
- Wide range of platforms for use with a variety of vessels.
- ♦ Maintenance- free brushless DC motor.
- Over speed detection and protection.



Specifications

Model	SK-L330-Pro
Cat. No.	8032120100
Amplitude	10mm
Max. load capacity (with platform)	7.5kg
Motor type	Brushless DC motor
Motor input power	28W
Motor output power	15W
Power	30W
Time settings range	1min-19h59min
Speed range	100-350rpm
Speed display	LCD
Timer display	LCD
Voltage	100-240V, 50/60Hz
Protection class	IP21
Dimension [W x D x H]	370×420×100mm
Weight	13.5kg
Permissible ambient temperature and humidity	5-40°C 80%RH
Data connector	RS232