

SK-0180-S

Features:

- Max.load capacity is 3kg with platform
- Double LED for independent display of time and speed
- Speed range of 40-200rpm
- Wide range of platforms for use with a variety of vessels



Specifications

Model	SK-O180-S
Cat. No.	8032210200
Orbital diameter	20mm
Max. load capacity (with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Power	20W
Time settings range	1min-19h59min
Speed range	40-200rpm
Speed display	LED
Timer display	LED
Voltage	100-240V,50/60Hz
Protection class	IP21
Dimension [W x D x H]	370×420×100mm
Weight	3.1kg
Permissible ambient temperature and humidity	5-40°C 80%RH