

Large Capacity Refrigerated centrifuge DL-6MC

Features:

- ♦ With the self-diagnostic and protection system for imbalance, over-speed, over-temperature. Greatly ensured the safe of users and machine.
- Microprocessor lid locks. The lid can be opened by electronic lock or via manual operation in case of resulted from broken electronic lock.
- Easy open the centrifuge lid via pneumatic springs;
- The centrifuge stand on movable casters.
- With the self-diagnosite system and protection for imbalance, over-speed and over-
- Temperature. Both automatic electronic lock and mechanical lock to ensure safety.

Application



Blood Bank



Bio-Technology



Chemistry Lab



Bio-Pharmaceutical



Parameters

Model	DI _6MC
Mar Cread	
Max Speed	6000r/min
Max RCF	6680xg
Max Capacity	6×1000ml
Speed Accuracy	±10r/min
Timer Range	0-99h 59 min
Temperature Range	-20~+40°C
Control&drive	Microprocessor control/direct drive
Acc/dec Rate	0-9(Grade)
Program Store	35
Motor Power	3.8KW
Refrigerator Power	2.2KW
Noise	≤60dB(A)
Power Supply	AC220V/110V 50Hz 35A
Weight(Without Rotor)	260kg
Dimensions	730×840×1200mm(L×W×H)
Refrigeration System	Big Brand High-performance Compressor CFC-free Refrigeration System
Rotor	No.1 AR 6*500ml (6000rpm 6680xg)
	No.2 SR 6*1000ml Round Cup (4200rom 5600xg)
	No.3 SR 6*1000ml Oval Cup (4000rom 5600xg)
	No.4 SR 222*5/7ml (4200rom 4390xg)