

## Large Capacity Refrigerated centrifuge DL-7MC

## Features:

- ♦ With the self-diagnostic and protection system for imbalance, over-speed, over-temperature. Greatly ensured the safe of users and machine.
- Microprocessor lid locks. The lid can be opened by electronic lock or via manual operation in case of resulted from broken electronic lock.
- Easy open the centrifuge lid via pneumatic springs;
- The centrifuge stand on movable casters.
- With the self-diagnosite system and protection for imbalance, over-speed and over-temperature.
- Both automatic electronic lock and mechanical lock to ensure safety.



## Parameters

Model	DL-7MC
Max Speed	8000r/min
Max RCF	11620xg
Max Capacity	6×1000ml
Speed Accuracy	±10r/min
Timer Range	0-99h 59 min
Temperature Range	-20~+40°C
Control&drive	Microprocessor control/direct drive
Acc/dec Rate	0-9(Grade)
Program Store	35
Motor Power	4KW
Refrigerator Power	2.2KW
Noise	≤60dB(A)
Power Supply	AC220V/110V 50Hz 35A
Weight(Without Rotor)	280kg
Dimensions	730×840×1200mm(L×W×H)
Refrigeration System	Big Brand High-performance Compressor CFC-free Refrigeration System
Rotor	No.1 AR 4*1000ml (8000rpm 11620xg) No.2 AR 6*1000ml (7000rom 11500xg)