

Lab Blenders Stomacher DW-09

Main Technical Parameters (DW-09)

1) Model: DW-09

2) Power: 300W

3) Supply: 220V/50Hz

4) Rattling time: 0.1~99 minutes (to 0.1minutes step) or 6~5940 seconds

5) Rap speed: 3~12 times/sec (to 0.1 times/second step)

Display the actual rattling number; the system automatically adjusts the slap speed

Multi-stage programming can be designed: parametric groups can be set according to one's wishes.

6) Parameter storage! More than 10 groups of multi-stage programming parameters are stored

7) Capacity: 3~400ml8) Display: digital display

9) Rattling Spacing: 0~5mm adjustable

10) Sterile bag: 170X300mm

11) Net weight: 18.5kg

12) Dimension: 420X245X340mm13) Packing size: 740X330X380mm

14) Soft-start mode: Gentle start

15) Pause Function: To ensure sample security, With auto-stop anti-trap feature Removable windows and adjustable slap board, toughened glass security door.

16) Slap Box: Stainless steel



Remark:

- 50 pcs germfree disposable homogenization bags are free of charge.
- Every machine with a stainless steel sample holder and three plastic clips(free of charge), as below picture:
- Optional accessories of LAB blenders
 - 1) Filter homogeneous bag
 - 2) Disposable sterile tube

