

DW-280B Electric Coal and Gas Heating Portable Autoclave

Features :

- High Quality SUS304 stainless steel material
- Electric, coal and gas heating
- Silicone seal
- Dual scale indicate pressure gauge
- Double safety valve to ensure safe
- Easy operation, safe and reliable



Parameters

Model	DW-280A	DW-280A	DW-280A
Volume	18L	24L	30L
Chamber size	Φ300*330mm	Φ300*460mm	Φ300*520mm
Working temperature	115-129°C	115-129°C	115-129°C
Working pressure	0.145-0.165Mpa	0.145-0.165Mpa	0.145-0.165Mpa
Packing size	40*40*45cm	40*40*53cm	40*40*64cm
G.W/N.W	16kg/15kg	17kg/16kg	18kg/117kg
Power	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW	AC220V/50Hz/2.0KW