


Horizontal Shaking Incubator


Features

- ◆ 1.Integrated with incubator and shaker to save space and cost.
- ♦ 2. High quality steel shell, polished stainless steel chamber.
- ♦ 3.Large LCD screen to display temperature and shaking speed.
- ◆ 4. With operation data memory function to eliminate faulty operation.
- ◆ 5.Non-volatile memory saves settings during a power outage and automatically restarts the unit as originally programmed after power is restored.
- 6.Auto-stop operation when door is opened. Strong air spring rod with easy opening and closing.
- ♦ 7.Brushless DC motor, more stable and reliable.
- 8. With UV Sterilizer
- 9. Equipped with leakage protection.









Accessory

Specification

Model	DW-SI-211B	DW-SI-211C
Shaking Speed	20-300rpm	
Speed Precision	±1 rpm	
Swing Amplitude	φ26mm	
Standard Configuration	500ml×28	2000ml×12
Max. Volume	250ml×36,500ml×28 ,1000ml×18	1000ml×18, 2000ml×12, 3000ml×8, 5000ml×6
Tray Size(mm)	920×510	920×510
Tray Included	1 piece	
Timing Range	1~9999min	
Display	LCD	
Temperature Range	RT+4~60°C(cooling)	
Temperature Accuracy	±0.1°C(constant temperature)	
Temperature Uniformity	±1°C	
Outer Size(cm)	120×74×80	120×74×100
Net Weight (kg)	174	183
Chamber Volume	970×565×280 (155L)	970×565×480 (265L)
Power Consumption	866W	951W
Power Supply	AC220V-240V 50~60HZ	