

Pharmacy Refrigerator

MPC-5V1006

Drawell 2~8°C pharmacy refrigerator are specially designed to store medicines,vaccines,regents and bio-medical products.

Features

- Accurate temperature controller, ensure the inside temperature within 2~8°C and accuracy 0.1 °C.
- Specially designed forced-air cooling system, ensure frost-free cabinet and good temperature uniformity.
- Suggested ambient temperature 10~32°C with humidity below 95%.
- Three layer toughened glass door with heater to avoid condensation.
- Color spray steel plate inner chamber for friendly using and easy cleaning.
- Wide available voltage range: 187V~242V

Safety

- Safety door lock and controller setting password to prevent unauthorized access.
- Build-in back-up battery to display temperature and to operate audible and visual alarm system for up to 8 hours without AC power.
- Auto on/off cooling fan to prevent cold air from leaking out when the door is opening.
- Visual and audible alarm system: High/low temperature, power failure, sensor error, door ajar
- Power failure protection: turn on delay of the cooling system after power failure.
- Controller error protection: the cooling system will remember the normal working cycle and keep working depending on this memory when controller or two sensor failed.

Humanized design

- 4 casters,and two with stopper for easy moving and fixation.
- 10 adjustable shelves for easy storing.
- Special designed service-free drainage system, ensure no water in the cabinet.
- Auto on/off LED light with manual switch on control panel.
- Double glass door automatic spring-back design

One year warranty, free accessories and technical support during the product warranty period



MPC-5V1006

Parameters

Performance

Temperature Range	2~8℃
Ambient Temperature	+10~+32℃
Controller	Microprocessor
Display	Digital Display

Refrigeration System

Compressor	1
Refrigerant	R600a
Cooling Type	Forced-air cooling
Defrost Type	Frost free

Alarm

High/low temperature	Y
Power failure	Y
Sensor error	Y
Door ajar	Y

Electrical Data

Power Supply(V/Hz)	220/50~60; 115/60
Rated Power(W)	472
Electrical Current(A)	2.71
Power Consumption	4.82kWh/24h

Dimensions

Capacity	1006L
Interior Size(W*D*H)	1100*684*1325mm
Exterior Size(W*D*H)	1220*872*1885mm
Package Size(W*D*H)	1330*950*2000mm
N.t./G.t.(approx)	189/213kg
Loading Qty(20'/40'/40'H)	8/16/16

Accessories

Caster	4
Lamp	LED
shelf	10
Temp. Recorder Chart recorder	Optional
Test hole	Y