

## Economic Vortex Mixer Model: MX-E

Vortex Mixers are ideal instruments for mixing liquid components in tubes, resuspension of cells, using an eccentric mechanism. They are widely used in biological and chemical analysis applications.

## Features:

- Touch operation
- Fixed speed at 3000rpm
- Small compact shape, ultra shock mount
- Eccentric bearing design



## Parameters

110-120V/220-230V
50/60Hz
Orbital
4.8mm
12W
10W
3000rpm
Inching
133×133×80mm
250×185×130mm
0.6kg
5-40°C
80%RH
IP21