

Tabletop PRP centrifuge TD4-PRP

TD4-PRP is Fat and PRP purified centrifuge, it can use 10ml, 20ml syringe to separate directly, the safe transfer process ensure the high survive rate of PRP and fat ,and reduced the complex operation process, it is the idea machine for cosmetic surgery, beauty salon, clinics and so on . can also fit many type of PRP kits in market.

Features:

- ◆ Achieved the maximum amount of PRP extraction through clinical research and test and special design on the Acc and Dec and RCF.
- Microprocessor control, Frequency conversion brushless AC motor.
- Touch LCD display with PRP program, easy to operate.
- Protection from imbalance, lid-lock, over-speed and display the error code on the screen.
- ♦ It is equipped with free time brake and able to set time, speed, Acc/Dec according to the requirement.
- It is equipped with free time brake and able to set time, speed, Acc/Dec according to the requirement.



Parameters

Model	TD4-PRP				
Max Speed/RCF	4000r/min / 3050xg				
NO.1 Swing Rotor	4*20ml syringe	4000r/min		3050xg	
NO.2 Swing Rotor	4*60ml PRO PRP tube	3158r/min		2810xg	
NO.3 Swing Rotor	Customize Made	4000r/min		3050xg	
Program	No.1	No.1		No.2-9	
	PRP			adjustable	
Motor	Microprocessor & Brushless DC Motor				
Timer Range	0~99min				
Noise	≤50dB(A)				
Power Supply	AC220V&110V 50Hz 2A				
Dimension	465*395*310mm(L*W*H)				
Weight	27 kg				