

UT710V Veterinary Infusion Pump

The UT710V Infusion Pump is an intelligent and precisely controllable infusion pump, which can be used in hospitals, clinical environments and emergency land transportation. 3000 mAh super long battery life, can be used continuously for 10 hours. Dual CPU and numeric keyboard control infusion is safer and more efficient. There are 5 modes to choose from to meet different needs.

Features:

- Dual CPU control infusion is safer.
- High-resolution LCD display.
- Numeric keyboard operation is more convenient and accurate.
- Dedicated blocking pressure sensor, dynamic display pressure value.
- 10 levels of blocking, bubble and alarm sound are adjustable.
- ◆ 3000mAh lithium battery backup life: ≥10hours.
- There is the first (pill) dose function, and the flow rate can be set arbitrarily.
- Start the night mode, after 5 minutes of no operation, the screen will automatically reduce the brightness.





High-resolution LCD screen, information is more intuitive.



The numeric keyboard makes the operation more convenient and precise.



The handle design confo- rms to ergonomics and is convenient to carry.

- Standard configuration: UT710V Main Unit/Power Cord
- Optional Configuration: Wireless Infusion Monitoring/Heater/Drop Sensor