

Thermostatic Floor Type Shaking Incubator

Applications

The product is suitable for medical and health care, pharmaceutical, biochemical, industrial production and agricultural science and scientific research departments for cell culture, fermentation, hybridization, enzyme and cell tissue research.

Features

- ◆ 1.Integrated with incubator and shaker to save space and cost.
- ◆ 2. High quality steel shell, polished stainless steel chamber.
- ◆ 3.Unique slope design panel with large LCD display, better visual effect.
- ◆ 4. With operation data memory and password locking function to eliminate faulty operation.
- ◆ 5.Non-volatile memory saves settings during a power outage and automatically restarts the unit as originally programmed after power is restored.
- ♦ 6. When the measured temperature deviates from set temperature by 38 (can be set
- ♦ 7.arbitrarily) will automatically stop heating and send out audible and visual alarms.
- ◆ 8.Equipped with strong and fast refrigeration system for instant cooling requirements, automatic defrost.
- ♦ 9. The most advanced large torque motor to ensure continuous work without maintenance.
- ◆ 10. With large view window and built-in lighting for easy observation.
- ♦ 11.Low noise, luxury and beautiful appearance.
- ◆ 12.Equipped with leakage protection.





DW-SI-2102 / DW-SI-2112B

Specification

Model	DW-SI-2102C(single)	DW-SI-2102(double)	DW-SI-2112B(double)
Rotation frequency (rpm)	40-300		
Frequency accuracy (rpm)	±1		
Shimmy amplitude(mm)	Ф26		
The standard configuration	250mlx24	250mlx35, 500mlx24	250mlx45, 500x37
The maximum capacity	100mlx48, 250mlx24, 500mlx18	250mlx70 ,500mlx48 1000mlx24, 2000ml×8	250mlx90, 500mlx74 1000mlx36, 5000mlx8
The tray number	2 pieces		
The pallet size(mm)	500x350	740x460	970x760
Power (W)	550	750	1420
Operating temperature	4~60°C	4~60°C	4~60°C
Temperature Accuracy	±0.1°C		
Temperature uniformity	±1°C		
Timing range	1-9999min		
Digital display mode	LCD		
The power supply	AC220V~240V±10% 50~60HZ		
Volume (cm) (W) x (D) x (H)	66x44x65 (170L)	84x52x64(280L)	109.5x62.5x85(582L)
Overall dimensions (cm)	70x56x126	95x70x126	120x82x160
Te net weight(KG)	165	240	425