



## D Series Ultrasonic Cleaning Machine

Digital display and setting ultrasonic cleaning time

Display working countdown time

Timing range: 1-600 minutes adjustable

Equipped with a special stainless steel basket, reduce sound cover

The shell and cover are made by stainless steel

The equipment controlled by Single-chip controller

The latest working parameters can be stored in memory after power-off

### Specifications:

Model	Inner Tank	Volume	Frequency	Power	Adjustable time	Drainage	Net weight	Packing size
	L/W/H ( mm)	L	KHz	W	min		Kg	mm
DW-50	150 × 85 × 65	1	40	50	1-120	/		220 × 150 × 210
DW-80	155 × 140 × 100	2	40	80	1-120	/	2.2	240 × 230 × 300
DW-100D	300 × 150 × 100	4.5	40	100	1-600	/	4.5	400 × 240 × 435
DW-120D	240 × 140 × 100	3	40	120	1-600	/	4.5	325 × 230 × 300
DW-120D	240 × 140 × 150	5	40	120	1-600	/	4	325 × 230 × 350
DW-3200D	300 × 155 × 150	6	40	180	1-600	/	5	400 × 240 × 380
DW-5200D	300 × 240 × 50	10	40	200	1-600	✓	6.8	450 × 325 × 375
DW-5200D	300 × 240 × 150	10	40	250	1-600	✓	7	450 × 325 × 375
DW-5200D	300 × 240 × 150	10	40	300	1-600	✓	7.2	450 × 325 × 375
DW4200D	300 × 240 × 200	14.4	40	400	1-600	✓	9	445 × 320 × 460
DW25-12D	500 × 300 × 150	22.5	40	500	1-600	✓	13	640 × 395 × 375
DW25-12D	500 × 300 × 150	22.5	40	600	1-600	✓	13.5	640 × 395 × 375

DW25-12D



DW-80



DW-50D



DW120D



DW-5200D

