



DTDN Series Ultrasonic Cleaning Machine

Large-screen LCD display

Display the parameters setting menu

Adjustable power: 40 %--100 %

Digital setting ultrasonic cleaning time

Working time countdown display

Display actual temperature of the cleaning tank

Temperature setting Range: room temperature--80°C

Timing Range: 1-600 minutes adjustable

Equipped with a special stainless steel basket, reduce sound cover

The shell and cover are made by stainless steel

The equipment controlled by Single-chip controller

The latest working parameters can be stored in memory after power-off

Model	Inner Tank		Volume	Frequency	power	Heating power	Acumulative time	Drainage
	L/W/H (mm)	L	KHz	W	W	min		
DW-120DTDP	240×140×100	3	40	120	150	99999	/	
DW-3200DTDN	300×155×150	6	40	180	150	99999	/	
DW-5200DTDN	300×240×150	10	40	200	800	99999	√	
DW-5200DTDN	300×240×150	10	40	250	800	99999	√	
DW-5200DTDN	300×240×150	10	40	300	800	99999	√	
DW-4200DTDN	300×240×200	14.4	40	400	800	99999	√	
DW25-12DTDN	500×300×150	22.5	40	500	1000	99999	√	
DW25-12DTDN	500×300×150	22.5	40	600	1000	99999	√	
DW-3200DTD	300×155×150	6	40	180	150	99999	/	
DW-5200DTD	300×240×150	10	40	200	800	99999	√	
DW-5200DTD	300×240×150	10	40	250	800	99999	√	
DW-5200DTD	300×240×150	10	40	300	800	99999	√	
DW-4200DTD	300×240×200	14.4	40	400	800	99999	√	
DW25-12DTD	500×300×150	22.5	40	500	1000	99999	√	
DW25-12DTD	500×300×150	22.5	40	600	1000	99999	√	

DW-3200DTDN



DW25--12DTDN

