



## DTS Series of Dual-frequency Ultrasonic Cleaning Machine

Interchangeable use of two kinds of frequencies: 25KHz, 40KHz

Digital display and setting ultrasonic cleaning time

Display working countdown time

Display actual temperature of the cleaning tank

Temperature setting range: Ambient temperature -- 80 °C

Timing range: 1-- 99 minutes adjustable

Equipped with a special stainless steel basket, reduce sound cover

The shell and cover are made by stainless steel

The equipment controlled by Single-chip controller

The lastest working parameters can be stored in memory after power-off

The equipment controlled by Single-chip controller

The lastest working parameters can be stored in memory after power-off

Model	Inner Tank	Volume	Frequency	power	heating power	Accumulative time	Drainage
	L/W/H ( mm)	L	KHz	W	W	min	
DW-3200DTS	300×155×150	6	25/40	180	120	1-99	/
DW-5200DTS	300×240×150	10	25/40	200	320	1-99	✓
DW-5200DTS	300×240×150	10	25/40	250	320	1-99	✓
DW-5200DTS	300×240×150	10	25/40	300	320	1-99	✓
DW-4200DTS	300×240×200	14.4	25/40	400	320	1-99	✓
DW25-12DTS	500×300×150	22.5	25/40	500	600	1-99	✓
DW25-12DTS	500×300×150	22.5	25/40	600	600	1-99	✓

DW-5200DTS



DW-4200DTS



DW25-12DT



DW-3200DTS

