

DW-B Digital Vertical Autoclave

Features :

- High quality SUS304 stainless steel material
- Dual scale indicate pressure gauge
- Light indicate working status
- LED digital display screen
- Microcomputer control system
- Temperature, time of sterilization can be set.
- Water lack auto-protection
- Over pressure auto-discharging
- Over current auto-protection
- Over temperature auto-protection
- Drying function can be add(best large than 75L)



Parameters

Model	DW-B35LDigital	DW-B50LDigital	DW-B75LDigital
Volume	35L	50L	75L
Chamber size	Φ300*485mm	Φ320*545mm	Φ388*480mm
Working temperature	100-134℃	100-134℃	100-134℃
Working pressure	0.22Mpa	0.22Mpa	0.22Mpa
Time scope	0-999min	0-999min	0-999min
Packing size	53*52*92cm	56*56*105cm	62*60*109cm
G.W/N.W	49kg/34kg	60kg/40kg	80kg/57kg
Power	AC220V/50Hz/2.5KW	AC220V/50Hz/2.5KW	AC220V/50Hz/3.5KW
Model	DW-B100L Digital	DW-B120L Digital	DW-B150L Digital
Volume	100L	120L	150L
Chamber size	Φ388*630mm	Φ460*590mm	Φ460*700mm
Working temperature	100-134℃	100-134℃	100-134℃
Working pressure	0.22Mpa	0.22Mpa	0.22Mpa
Time scope	0-999min	0-999min	0-999min
Packing size	63*61*118cm	66*67*121cm	67*68*128cm
G.W/N.W	85kg/66kg	97kg/75kg	110kg/80kg
Power	AC220V/50Hz/4.0KW	AC220V/50Hz/4.0KW	AC220V/50Hz/4.0KW