

0.01g YP Electronic Balance LED Display (External Calibration)



• High precision loadcell.

- Rugged Plastic Housing.
- Big size Stainless Steel Platter.
- High contrast easy to read LED display.
- Built-in rechargeable battery last up to 7 days.
- Height adjustable feet.
- Full capacity substraction.
- Stabilization time 1.5 seconds.
- Overload Protection.
- Below Balance Weighing Facility.
- Checking Weighing, piece counting function.
- Can be contain bubble level.



Specifications

Model	YP1002	YP2002	YP3002	YP5002	YP6002	YP10002	YP20002
Capacity	0-100g	0-200g	0-300g	0-500g	0-600g	0-1000g	0-2000g
Readability	0.01g						
Repeatability	±0.02g						
Linearity	±0.02g						
Platter Size	143*192mm/145*205mm						
Dimensions	280mm x 180mm x80mm						
Net Weight	1500g						



www.drawell.com.cn