DRAWELL Artist of Science

0.1g /1g YP Large Balance LED



- High precision loadcell.
- Rugged Plastic Housing.
- Big size Stainless Steel Platter.
- High contrast easy to read LED display.
- Built-in rechargeable battery last up to 7 days.
- Height adjustable feet.
- Full capacity substraction
- Stabilization time 1.5 seconds
- Panel simple para una facil operation
- Overload Protection
- Below Balance Weighing Facility
- Checking Weighing, piece counting function.
- RS232 Interface.



Specifications

Model	YP60001L	YP100001L	YP150001L	YP200001L	YP300001L	YP15000L	YP30000L	
Capacity	0-6000g	0-10000g	0-15000g	0-20000g	0-30000g	0-15000g	0-30000g	
Readability	0.1g					1g		
Repeatbility	±0.2g	±0.2g					±2g	
Linearity	±0.2g		±2g					
Pan size	335*235mm							
Dimensions	370*370*125	Smm						
Net Weight	5000g							

	Power	Supp	ļ
--	-------	------	---

www.drawell.com.cn