

0.1g YP Electronic Balance (External Calibration)



• High precision loadcell.

- Rugged Plastic Housing.
- Big size Stainless Steel Platter.
- High contrast easy to read LED display.
- Built-in rechargeable battery last up to 7 days.
- Height adjustable feet.
- Full capacity substraction
- Stabilization time 1.5 seconds
- Panel simple para una facil operacion
- Overload Protection
- Below Balance Weighing Facility
- Checking Weighing, piece counting function.





Specifications

Model	YP2001D	YP6001D	YP10001D	YP15001D	YP20001D	YP30001D	YP50001D
Capacity	0-200g	0-600g	0-1000g	0-1500g	0-2000g	0-3000g	0-5000g
Readability	0.1g						
Repeatbility	±0.2g						
Linearity	±0.2g						
Pan size	335*235mm						
Dimensions	370*370*125	Smm					
Net Weight	5000g						

	Power	Supp	ļ
--	-------	------	---

www.drawell.com.cn