

## DW-LDO-A Series Forced Air Drying Oven

## Reatures

- Vertical forced air type good temperature uniformity
- Microprocessor controller (with temperature correction and timing function), accurate and reliable.
- Large LCD screen display.
- High quality stainless steel chamber, removable shelf, easy-to-clean.
- Silicon sealing ring for reliable sealing.
- Equipped with leakage protection.
- Equipped with spare temperature control Ensuries to work normally when the
- main temp control failed.
- Anti-hot handle.



● Temp. Fluctuation:±1°C.

• Temp. Uniformity:±2.5°C%.

Temp. Resolution: 0.1°C

Timer: 1~9999m

## Optional Parts

- Multi-segment programmable control
- Built-in Printer
- RS485 interface



## Specifications

Model	Volume(L)	Chamber Size (WxDxH)cm	Power Rating (W)	Temp. Range (°C)	Shelf/Max	Net/Gross(kg)	Power Supply
DW-LDO-9023A	25	30×30×28	320	RT+5~200	2/4	29/36	AC 220V 50HZ
DW-LDO-9053A	54	42x37x35	600		2/5	35/43	
DW-LDO-9123A	138	50x50x55	120		3/7	45/55	
DW-LDO-9203A	252	60×60×70	1600		3/10	60/69	
DW-LDO-9030A	31.5	30×30×35	700		2/5	29/37	
DW-LDO-9070A	72	40×40×45	1000		2/7	45/52	
DW-LDO-9140A	136	45×55×55	1600		3/8	45/56	
DW-LDO-9240A	225	50×60×75	2000		3/11	62/70	
DW-LDO-9420A	429	60×55×130	2600		4/6	130/153	AC 380V 50HZ
DW-LDO-9620A	624	80×60×130	3600		5/6	168/186	
DW-LDO-9036A	31.5	30×30×35	950	RT+5~300	2/5	29/37	AC 220V 50HZ
DW-LDO-9076A	72	40×40×45	1400		2/7	40/50	
DW-LDO-9146A	136	45×55×55	2000		3/8	45/56	
DW-LDO-9246A	225	50×60×75	2600		3/11	62/70	
DW-LDO-9426A	429	60×55×130	3500		4/6	130/153	AC380V 50HZ
DW-LDO-9626A	624	80×60×130	5700		5/6	168/186	