

CGF-iPRF-APRF-PRP centrifuge TD4A

TD4A is PRP (platelet-rich Plasma), PRF (platelet-rich fibrin) and (CGF)Concentrate Growth Factors centrifuge, It's widely used in Dental hospital, Plastic surgery, Orthopedics department, and Maxillo facial surgery.

Features

- 1. With compact structure, low noise and light weight
- 2. Large LCD display, with touch panel. Easy to operate
- 3. High efficiency for extract.
- 4. With dedicated program for PRP, iPRF, APRF, CGF, also can used for ordinary medical use
















Parameters

Model	TD4A						
Max Speed		4000r/min					
Rotor	Capacity		Max Speed		M	Max RCF	
NO.1 Angle Rotor	6*10/15ml		4000r/min		18	1880xg	
NO.2 Angle Rotor	8*10/15ml			4000r/min	18	1800xg	
Program	No.1	No.2		No.3	No.4	No.5	
	CGF	PRP		A-PRF	iPRF	adjustable	
Motor	Microprocessor & Brushless DC Motor						
Timer Range	0~99min						
Noise	\leq 50dB(A)						
Power Supply	AC220V&110V 50Hz 2A						
Dimension		315*270*215mm(L*W*H)					
Weight		5kg					