

## Benchtop Automatic Balance Low-Speed centrifuge TDZ5-WS

TDZ5-WS low speed centrifuge are widely used in medical and Biological Field, Health Examination center, and laboratories etc.

## Features:

- Stainless Steel Chamber, Full metal structure and microprocessor controlled
- ◆ Safety Interlock system that prevents opening the centrifuge lid during running.
- With Emergency lid lock release, in case of power off.
- Imbalance sensor, Over speed
- 10 programs can be stored, quickly recall the last used program.
- With air cooling fans to avoid high temperature with long time working.
- Digital or LCD display for choose, which indicated the Speed, time, RCF etc.





Model	TDZ5-WS		
Max Speed	5500r/min		
Speed Accuracy	$\pm 10 rpm$		
Max RCF	5220×g		
Timer Range	0~99mins		
Noise	≤65dB(A)		
Power Supply	AC220V & 110V 50Hz 5A		
Dimension	$580\times450\times360$ mm(L×W×H)		
Weight	39kg		
Rotor No.	Rotor Capacity	Speed (r/min)	RCF (xg)
No.1 Angle Rotor	12×10/15ml	5500	5220
Vo.2 Swing Rotor	4×100ml	4000	3020
No.3 Swing Rotor	8×50ml	4000	3020
	8×100ml	4000	3020
No.4 Swing Rotor	32×10/15ml	4000	3020
No.5 Swing Rotor	24×10/15ml	4000	3020
No.6 Swing Rotor	16×10/15ml	4000	3020
No.7 Swing Rotor	48×5ml	4000	3020
No.8 Swing Rotor	4×250ml	4000	3300
	48×10/15ml	4000	3000
No.9 Swing Rotor	72×2/5ml		2300