

## DW-SXT-02 Soxhlet Extractor

DW-SXT-02 Soxhlet Extractor is for labs to make Soxhlet extraction, designed in accordance with the classical Soxhlet extraction method. This apparatus is reasonably designed, provides stable performance, high precision, time-saving operation, is an ideal Soxhlet extraction apparatus for foodstuff, oil&fat, feedstuff and other industries.

### Features

- A. Take water bath to heat, providing more even warming;
- B. Solvent for extraction can be recovered and used repeatedly;
- C. Extraction duration can be set up at random, alarm when time is out.

### Technical Parameters

- 1. Measurement range: for extract samples like cereals, feedstuff, oily seed, soil, etc.
- 2. Sample number:2/batch
- 3. Capacity of extraction flask:500ml/per
- 4. Sample amount:0.5-20g/per
- 5. Temp. range:ambient+5°C~100°C
- 6. Power: 220V+10V, 50Hz
- 7. Heating power:300W
- 8. Dimensions(mm):750×360×550
- 9. Weight:16kg

