

DW-LDO Constant Temperature Drying oven

Features

High quality stainless steel chamber, safe and easy-to-clean.

Microprocessor controller with LCD display, accurate and reliable.

Equipped with multiple sets of heaters.

- Equipped with leakage protection.
- Equipped with spare temperature control which ensures the product work normally even the main temp.control failed.
- Adopt temperature selection switch which can select the grade according to the heating speed and operating temperature.
- Convenient operation for air inlet, reasonable air duct structure, good temperature uniformity.
- With inner glass door for easy observation.
- Anti-hot handle

Technical Parameters

- Temp. Range:RT+10°C-300°C
- Temp. Fluetuation: $\pm 1^{\circ}$ C.
- Temp.Resolution: 0.1°C
- Temp. Uniformity: \leq Max. Temp. $\pm 3.5^{\circ}$ C%.



Optional Parts

- Multi-segment programmable control
- Built-in Printer
- RS485 interface

Specifications

Model	Volume(L)	Chamber Size (LxWxH)cm	Power Rating (W)	Shelf	N.W/G.W (kg)
DW-LDO-300	43	35×35×35	1300	2	35 / 55
DW-LDO-400	81	45x40×45	2000	2	50 / 80
DW-LDO-500	138	50×50×55	2500	2	75 / 110
DW-LDO-600	252	60×60×70	3500	2	110 / 150

