

JY88-IINs JY96-IIN Ultrasonic Homogenizer

JY88-IIN、 JY96-TIN Ultrasonic Homogenizer, also known as ultrasonic emulsification machine, ultrasonic dispersion machine, should not be dealing with micro-sample



Parameters

MODEL	JY88-IIN	JY96-IIN
Frequency	20-25 KHz	20-25 KHz
Power	250 W (1%- 100%)	150 W (1%- 100%)
Accessories (Φ)	Φ 6mm	Φ6mm
Optional Tip (Φ)	Φ 3mm	Φ 3mm
Crushing Capacity	10-200 ml	10- 100 ml
Temperature alarm range	0-99.9C	0-99.9C
Alarm items	Time, temperature, overload	Time, temperature, overload
Duty Ratio (%)	1-99 %	1-99 %
Supply	220/110V 50Hz/60Hz	220/110V 50Hz/60Hz
Power supply chassis Size and Weight	430*245*300mm; 8.5Kg	430*245*300mm; 8.4Kg
Sound chamber size&weight	345*345*535mm; 12kg	345*345*535mm; 12kg