

JY99-IIDN Ultrasonic Homogenizer



Parameters

| MODEL No | JY99-IIDN |
|----------------------------|-----------------------------|
| Frequency | 19.5-20.5KHz |
| Power | 1800W (1%- 100%) |
| Display | 7" TFT Touch display |
| Temperature alarm range | 0-99.9C |
| Accessories (Φ) | Φ25 |
| Working method | Interval /continue |
| Alarm items | Time, temperature, overload |
| Optional Tip (Φ) | Φ15mm, 28mm |
| Crushing Capacity | 50-1200ml |
| Duty Ratio (%) | 0.1-99.9% |
| Supply | 220/110V 50Hz/60Hz |
| Host size& wight | 430mm*245mm*300mm &11.8Kg |
| Sound chamber size &weight | 345mm*345mm*535mm &12kg |