

SK-D1807-S

Features:

- ♦ Three-dimensional rocking motion with the angle of 7°
- Max. load capacity is 3kg with platform
- Speed range of 10-80rpm
- Double LEDs for independent display of time and speed
- Various of platform accessories



Specifications

Model	SK-D1807-S
Cat. No.	8032412000
Tilt angle	7°
Max. load capacity (with platform)	3kg
Motor type	DC motor
Motor input power	16W
Motor output power	10W
Speed range	10-80rpm
Speed display	LED
Timer display	LED
Time settings range	1min-19h59min
Voltage	100-240V, 50/60Hz
Power	20W
Weight	2.7kg
Permissible ambient temperature and humidity	5-40°C 80%RH
Dimension[W×D×H]	270×330×130mm
Protection class	IP21