

DWS-280 Series Portable Pressure Steam Sterilizer

Features

A, A+ type Electric heating type

C,C+type, Electric or LPG heating type

1. High quality stainless steel material
2. Simple operation, safe and reliable
3. Electric heating mode
4. Auto release when over pressure
5. Double scale indication pressure gauge

B,B+ Self control electric heating type

1. All made of SUS304 stainless steel material, durable
2. Microcomputer control,digital display,sterilization temperature and time can be adjusted
3. Auto release when over pressure,Low water protection,prevent dry burning
4. Double scale indication pressure gauge, clear to read
5. Simple operation,safe and reliable



DWS-280A/A+



DWS-280B/B+



DWS-280C/C+

Parameters

Mode	DWS-280A	DWS-280A+	DWS-280B	DWS-280B+	DWS-280 C	DWS-280 C+
Volume	18L	24L	18L	24L	18L	24L
Power	2.0kw	2.0kw	2.0kw	2.0kw	2.0kw	2.0kw
Voltage	110V/220V, 60Hz/50Hz				AC 220V、 50Hz	
Rated working pressure	0.142MPa					
Rated working temperature	126°C					
Chamber volume	Ø 280×295	Ø280×390	Ø280×295	Ø280×390	Ø280×295	Ø280×390
Packing size(LxWxH)	410×410×560	410×410×680	410×410×560	410×410×680	410×410×430	410×410×540
Gross/net weight	17/16kg	19/17.5kg	18/17kg	19/17.5kg	15.5/14.5kg	7.5/15.8kg