

## DWT-A Series Table Top Steam Sterilizer

### Features

1. Quality stainless steel inner tank, sealed cover
2. The microcomputer controls the automatic cycle, 121 °C and 134 °C can be set
3. Touch buttons
4. Equipped with standard test interface
5. European CLASS N standard.
6. With automatic protection function: over temperature protection, over pressure protection, low water level protection, anti dry burning
7. Buzzer reminder after sterilization, automatic stop.
8. Automatic discharge of cold air, automatic steam exhaust after sterilization.
9. Built-in steam generator to provide saturated steam quickly



### Parameters

Mode	DWT-18A	DWT-23A
Volume	18L	23L
Power	1.5kw	2kw
Voltage	110V/220V, 60Hz/50Hz	
Design pressure	0.25MPa	
Design temperature	139°C	
Rated working pressure	0.22MPa	
Rated working temperature	134°C	
Chamber volume	Ø250×370	Ø250×452
Outer dimension	500×400×350	600×400×350
Packing size(LxWxH)	540×450×390	640×450×390
Gross/net weight	26/24kg	29/27kg