

## EC10 Portable Cond Meter

### Features

- 3.7" Color screen, the user-friendly interface.
- Data logging up to 3000 groups.
- Calibration points: 1-pt.
- Equipped with a USB Type-C charging port, conforms with the EU new law.
- Ideal for use in harsh environments — IP57 Waterproof rating, silicone socket protector and grip cover protect the instrument from general impact.
- Large backlit LCD displays conductivity/TDS/salinity and temp along with icons of stable readings.
- Foldable stand allows use as a benchtop meter.



### Specifications

Cond	Range	0-2000.0mS/cm (auto divided by range)
	Accuracy	±0.5%FS ± one digit
	ATC range	( 0 - 50)°C
Temperature	Range	0-100°C
	Resolution/Accuracy	0.1°C; ±0.5°C±1digit
Data storage/content	3000 groups (Serial number, Mode, Measuring value, Temp, Date: Y/M/D/H/M)	
Data export/Charging port	Type-C - USB	
PC software function	GLP measurement management Two-level user management The user's operation history log Meter usage management Automatically /manually timing data storage	
Power supply	Rechargeable Lithium battery × 1	
Packing including	1. Meter× 1 2. Plastic Cond electrode Con201T× 1 (K=1, PC housing, for general liquid test) 3. Standard Cond solution: 84uS/1413uS/12.88mS/cm, 50mL, one for each 4. USB-C Charger Cable× 1 5. Hard carrying case×1	