

TD5-WS


Intelligent Automatic Balance Low-Speed Centrifuge

Introduction

TD5-WS intelligent low speed centrifuge are widely used in medical and Biological Field, Flealth Examination center, and laboratories etc.

Features:

1. Stainless Steel Chamber, Full metal structure and microprocessor controlled
2. With electric door lid open and closed system, also equipped with emergency lid lock release in case of power off unexpected.
3. With the self-diagnostic and protection system for imbalance, over-speed, over-temperature.
Greatly ensured the safe of users and machine.
4. 10 programs can be stored, quickly recall the last used program.
5. With air cooling fans to avoid high temperature with long time working.
6. Full colored touch panel, shows with RPM, RCF, Acc, Dec , Rotor No., etc.
7. Automatic rotor identification avoid over speed running.



Parameters

Model	TD5-WS		
Max Speed	5500r/min		
Speed Accuracy	±10rpm		
Max RCF	5220xg		
Timer Range	0~99min		
Rotor No.	Rotor Capacity	Speed (r/min)	RCF (xg)
No.1 Angle Rotor	12x10/15ml	5500	5220
No.2 Swing Rotor	4x100ml	4000	3020
No.3 Swing Rotor	8x50ml	4000	3020
	8x100ml	4000	3020
No.4 Swing Rotor	32x10/15ml	4000	3020
No.5 Swing Rotor	24x10/15ml	4000	3020
No.6 Swing Rotor	48x5ml	4000	3020
No.7 Swing Rotor	72x2/5ml	4000	3020
No.8 Swing Rotor	4x250ml	4000	3300
No.9 Swing Rotor	48x10/15ml	4000	3000
No.10 Micro plates	2x2x96holes	4000	2300
Timer Range	0~99mins		
Noise	≤55dB(A)		
Power Supply	AC220V & 110V 50Hz 5A		
Dimension	600x460x340mm(LxWxH)		
Weight	39kg		