

TD6-WS

Intelligent Benchtop Large Capacity Low-Speed Centrifuge

Introduction

TD6-WS is low speed Large Capacity centrifuge, Which is the idea machine for Biological, hospital clinics laboratories, Epidemic prevention stations and to separate samples of blood, plasma, Urine and so on.

Features:

1. Stainless Steel Chamber, Full metal structure and microprocessor controlled
2. With electric door lid open and closed system, also equipped with emergency lid lock release in case of power off unexpected.
3. With the self-diagnostic and protection system for imbalance, over-speed, over-temperature. Greatly ensured the safe of users and machine.
4. 10 programs can be stored, quickly recall the last used program.
5. With air cooling fans to avoid high temperature with long time working.
6. Full colored touch panel, shows with RPM, RCF, Acc, Dec, Rotor No., etc.
7. Automatic rotor identification avoid over speed running.



Parameters

Model	TD6-WS		
Max Speed	6000r/min		
Speed Accuracy	±30rpm		
Max RCF	5600xg		
Timer Range	0~99min		
Rotor No.	Rotor Capacity	Speed (r/min)	RCF (xg)
No.1 Angle Rotor	12×10/15ml	6000	5600
No.2 Swing Rotor	24×5ml	4200	3120
No.3 Swing Rotor	108×5ml	4000	4200
No.4 Swing Rotor	16×10/15ml	5000	4100
No.5 Swing Rotor	8×50ml	4200	3120
	4×50ml		
No.6 Swing Rotor	8×100ml	4200	3120
No.7 Swing Rotor	4×500ml	4000	3600
No.8 Well Plates Rotor	2×2×96holes	4000	2300
Timer Range	0~99mins		
Noise	≤65dB(A)		
Power Supply	AC220V & 110V 50Hz 5A		
Dimension	605×505×420mm(LxWxH)		
Weight	52kg		