


TDL-6MX

Intelligent Benchtop Low Speed Refrigerated Centrifuge

Features:











1. Microprocessor control, speed control by frequency-conversion and CFC compressor. Maintenance free Frequency-conversion motor without carbon dust pollution.
2. Advanced door lid control system can both electric control open and close door lid, also equipped with emergency lid lock release to avoid power off unexpected.
3. 10 programs can be stored, quickly recall the last used program.
4. Automatic rotor identification to avoid over speed running by mistake operate.
5. Easy operation, 7 inches LCD touch panel, shows all parameters in screen.
6. With Emergency lid lock release, in case of power off.
7. Error minder for wrong operate or problem.



Parameters

Model	TDL-6MX	Speed Accuracy	±50r/min
Max Speed	6000r/min	Timer Range	0~99min
Max RCF	5600xg	Temp Range	-20°C~ +40°C
Display	LCD Touch display	Max Capacity	4x250ml
Noise	≤60dB(A)	Power Supply	AC220v&110V 50Hz 18A
Weight	85kg	Dimension	720x500x420mm

Rotor

 No.1 AR	Max Speed:6000r/min Capacity:12x15/10ml Max RCF:5600xg	 No.2 SR	Max Speed:4000r/min Capacity:24x10/15ml Max RCF:3020xg
 No.3 SR	Max Speed:4000r/min Capacity:32x10/15ml Max RCF:3020xg	 No.4 SR	Max Speed:4000r/min Capacity:4/8x50ml Max RCF:3020xg
 No.5 SR	Max Speed:4000r/min Capacity:72x2/5ml Max RCF:3020xg	 No.6 SR	Max Speed:4000r/min Capacity:48x5ml Max RCF:3000xg
 No.7 SR	Max Speed:4000r/min Capacity:4/8x100ml Max RCF:3020xg	 No.8 SR	Max Speed:4000r/min Capacity:4x250ml Max RCF:3300xg
 No.9 SR	Max Speed:4000r/min Capacity:48x10/15ml Max RCF:3020xg	 No.10 MTP	Max Speed:3000r/min Capacity:2x2x96holes Max RCF:2300xg