

# Tower Electric Heating Stainless Steel Distiller

## TZ50/TZ100/TZ200/TZ400/TZ600

### Performance and characteristic:

- All of products are made from stainless steel.
- Use high temperature steam to heat, save the energy.
- Use the steam condensation water.
- Seamless stainless-steel pipes, heat efficiency is better.
- Condenser is made of seamless stainless-steel pipes and so with very good cooling effects & great capacity of water produce.
- Percolation, pushing ammonia, dividing water and steam when distilling, in order to control the quality.



### Parameters

Model	Water discharge	Pressure of steam supply	Pressure of cooling water	Steam consumption	Ratio of water consumption	Product size (Wooden)	Gross Weight (kg)	Net Weight (kg)
TZ50	≥50L/h	0.4MPa	0.2MPa	≈1:1	≤1:14	75 X65 X145	88	55
TZ100	≥100L/h					77 X 67X 165	100	65
TZ200	≥200L/h					102 X 82X140	150	95
TZ400	≥400L/h					113 X103 X165	200	130
TZ600	≥600L/h							