

## Integrated All Steel Fume Hood

**Model: DW-FH1200 / DW-FH1500 / DW-FH1800**

The fume hood is a common laboratory equipment used in the laboratory, which needs to discharge harmful gases, and needs to be cleaned and drained during the experiment.

### Features

The water, electricity, gas and ventilation system are integrated, and the multi-function power socket is built in, which is convenient for the use of other electrical equipment in the laboratory process. The quick-opening valve is used to facilitate the convenience of water during the experiment. The front baffle is a glass door that can move up and down, and the top is a low-speed exhaust fan, which can smoothly discharge harmful and odorous gases during the experiment. The bottom of the working surface is equipped with a stainless steel water tank (optional), which can discharge the disinfectant and experimental residues from a drainage tank through water washing, so as to keep the experimental environment safe and reliable.

### Structure

The structure is seamlessly welded with cold-rolled steel plates, the tabletop in the experimental area is veneered with stainless steel, and the surface of the equipment is electrostatically sprayed. The appearance is smooth and dust-free, beautiful and durable.



### Parameters

型号 Model	DW-FH1200	DW-FH1500	DW-FH1800
参数 Spec.			
排风速度 Exhaust speed	0.25 ~ 0.45m/s范围可调, 0.25-0.45m/s, adjustable		
进水压力 Inlet pressure	> 0.5Pa		
电源 Power	AC单相220V/50Hz		
工作区尺寸 Work area size, W1×D1×H1	1030*660*680	1300*660*680	1600*660*680
外型尺寸 Overall size, W×D×H	1200*740*1950	1460*740*1950	1760*740*1950
LED lamp specifications & qty	28W*	30W*	40W*