

DEB Electronic Analytical Balance 1mg

Advantages

1. Electromagnetic force sensor
2. Super bright detachable displayer
3. Multi-unit weighing function, customizable units
4. RS232 and USB bidirectional interface
5. Glass windshield-antistatic material
6. Fast sensitivity setting
7. Below weighing system
8. Time display
9. The weighing results can be exported to Excel
10. Built-in auto calibration system(Intemal calibration model)

Weighing function

1. Basic weighing
2. Dynamic weighing
3. Piece counting
4. Percent weighing
5. Density determination
6. Pipette calibration function



Parameters

EX-Cal Model	DEB503	DEB603	DEB1003	DEB2003
IN-Cal Model	DEB503N	DEB603N	DEB1003N	
Capacity	0-520g	0-620g	0-1000g	0-2000g
Readability	1mg			
Stable time	2s			
preheating	30-6mins			
Weighing units	g/oz/ct/mg			
Min Weight	1mg			
Repeatability	±1mg			
Linearity	±2mg			
Interface	RS232 / USB			
Pan size	Φ90mm			
Packing size	Inner carton:430x320x415mm, Master carton:490x380x485mm			
N.W.	6280g			
Power	AC110V-240V			