


DHP-9015A/B

Electric Constant-Temperature Incubator

Features

- The product shell is made of high-quality cold plates and processed into shape. The surface is treated by spraying process.
- The working chamber can be manufactured from stainless steel.
- The intelligent digital PID temperature control instrument is equipped.
- The large double-layer tempered glass observation window is available, facilitating the observation of changes in the items inside the box.
- A new type of synthetic silicone sealing strip is adopted, which is aesthetically pleasing and can operate at high temperatures for a long time, with a longer service life.



Specifications

Model	DHP-9015A	DHP-9015B
Capacity	15L	15L
Power supply	220V 50Hz	220V 50Hz
Temperature control range	RT+5-65°C	RT+5-65°C
Temperature resolution	0.1°C	0.1°C
Power	180W	180W
Studio size (in mm)	250*250*250	250*250*250
Overall size (in mm)	355*470*400	355*470*400
Packaging size (in mm)	460*530*490	460*530*490
Gross weight	21kg	21kg
Inner liner material	Galvanized inner liner	201 stainless steel inner liner