

## DEB Electronic Analytical Balance 0.1mg

### Advantages

1. Electromagnetic force sensor
2. Super bright detachable displayer
3. Multi-unit weighing function, customizable units
4. RS232 and USB bidirectional interface
5. Glass windshield-antistatic material
6. Fast sensitivity setting
7. Below weighing system
8. Time display
9. The weighing results can be exported to Excel
10. Built-in auto calibration system(Intemal calibration model)

### Weighing function

1. Basic weighing
2. Dynamic weighing
3. Piece counting
4. Percent weighing
5. Density determination
6. Pipette calibration function



### Parameters

EX-Cal Model	DEB124	DEB224	DEB324	DEB424	DEB524
IN-Cal Model	DEB124N	DEB224N	DEB324N	DEB424N	DEB524N
Capacity	0-120g	0-220g	0-320g	0-420g	0-520g
Readability	0.1mg				
Stable time	1.5s		2s		
preheating	30-60mins				
Weighing units	g/oz/ct/mg				
Min Weight	0.1mg				
Repeatability	±0.1mg		±0.2mg		
Linearity	±0.2mg		±0.3mg		
Interface	RS232 / USB				
Pan size	Φ90mm				
Packing size	Inner carton:430x320x415mm, Master carton:490x380x485mm				
N.W.	6280g				
Power	AC110V-240V				