

## KFSH-01 Karl Fischer Oven Sampler



### GENERAL FEATURES

- Color High contrast LCD touchscreen, 7.0 inches.
- Support setting temperature, flow rate and min extraction time.
- It shows the instrument status, operation status (temperature, flow, etc.)
- Support indicator prompt, buzzer prompt.

Model		KFSH-01
Flow Rate	Range	(0~200)mL/min
	Resolution	1mL/min
	Accuracy	±10mL/min
Temperature	Range	(40.0~300.0)°C
	Resolution	0.1°C
	Accuracy	±5°C
Work environment		Ambient temperature: (0~40)°C Relative humidity: not more than 85%
Power Supply		AC(220±22)V, (50±1)Hz
Dimensions		360×240×350mm
Net Weight		About 4.5kg