

KFT-40C Karl Fischer Titrator



FEATURES

- LCD display screen.
- Support KF Titration Mode and Titer Detection Mode.
- Features of auto-filling, auto-purging and auto-mixing of the reagents ensure safe handling of Karl Fischer chemicals.
- Selectable units including mg, mg/L, %, ppm, etc.
- Settable parameters, including measurement unit, polarization current, stirring rate, titration rate, stop volume, endpoint potential, stop criterion and etc.
- Store up to 200 titration data sets (GLP-compliant).
- Data can be easily transferred to printer via RS-232 communication interface.
- Reset feature automatically resumes all settings back to factory default options.

	Model	KFT-40C
Volumetric Titration	Range	0.1mg~250mg
	Resolution	0.1mg
	Working Current	1 μ A \pm 0.2 μ A; 50 μ A \pm 10 μ A
	Repeatability	\leq 0.5%
General	Power Supply	AC Adapter, 100-240 V AC input
	Dimensions	340 \times 400 \times 400 mm
	Weight	10kg