

Ti-20 Potential Titrator



FEATURES

- LCD display screen.
- Support titration modes: SET (Preset Endpoint Titration) and MAT (Manual titration).
- Glass burette can be used for volume measurement (prepared by user).

Model		Ti-20
Measuring Unit	Range	(0~±1400) mV, (0~14.00)pH
	Resolutions	1mV, 0.01pH
	Accuracy	pH: ±0.03pH, mV: ±0.35%FS
	Fluctuation	±0.01pH /3h
General	Power Supply	AC Adapter, 200-240 VAC input
	Dimensions	300×235×100 mm
	Weight	5 kg