

Packing List

KFT-41V Karl Fischer Titrator

No.	Description	Quantity
1	KFT-41V Karl Fischer Titrator	1 pc
2	Solvent manager	1 pc
3	500mL solvent bottle including cap and GL45 connector	1 pc
4	Desiccator tube for solvent bottle	1 pc
5	1000mL waste bottle including cap and GL45 connector	1 pc
6	500mL titrant bottle including cap and GL45 connector	1 pc
7	Volumetric vessel(including septum, seal, stopper and screw cap)	1 pc
8	CDY-4 Measuring electrode	1 pcs
9	Stirring bars($\varphi 7 \times 30$)	3 pcs
10	Drying cylinder	1 pc
11	Air tubing for waste(L=520mm, including joints)	1 pc
12	Air tubing for fill(L=580mm, including joints)	1 pc
13	Drain tube (L=800mm, including joints)	1 pc
14	Intake tube (L=800mm, including joints)	1 pc
15	Titration inlet tube (L=700mm, including joints)	1 pc
16	Pump tube (L=150mm, including joints)	1 pc
17	Cylinder for oven sampling tube	1 pc
18	Micro injector(100 μ L)	1 pc
19	Micro injector(1mL)	1 pc
20	Septum($\varphi 13 \times 3$)	10 pcs
21	Ejector tool (SC8.227.126)	1 pc
22	Titration outlet tube(L=450mm, including joints) and titration tip with anti-diffusion valve	1 pc
23	Burette device 10mL	1 pc
24	Burette cover	1 pc
25	Fuse $\Phi 5 \times 20$ (4A)	3 pcs
26	Double magnetic ring power cord (black)	1 pc
27	M12-5 cable(Connect the Karl Fisher Titrator to the solvent manager)	1 pc

No.	Description	Quantity
28	RS232-DB9 9 Pin male to male Serial Cross Cable-1 Meter	1 pc
29	USB communication cable-1.8 Meter	1 pc
30	USB to RS232 DB9 Serial Port Cable-1.5 Meter	1 pc
31	Touch pen (Capacitive screen)	1 pc
32	Instruction manual	1 pc
33	Sealing grease	1 pc
34	5A molecular sieve(500g)	1 pc