

Packing List

KFT-50VC Karl Fischer Titrator

No.	Description	Quantity
1	KFT-50VC Karl Fischer Titrator	1 pc
2	Solvent manager	1 pc
3	500mL solvent bottle including cap and GL45 connector	1 pc
4	Desiccator tube for solvent bottle	1 pc
5	1000mL waste bottle including cap and GL45 connector	1 pc
6	500mL titrant bottle including cap and GL45 connector	1 pc
7	Volumetric vessel (including septum, seal, stopper and screw cap)	1 pc
8	Coulometric vessel(including septum, stopper and screw cap)	1 pc
9	CDY-4 Measuring electrode	2 pcs
10	CDJ-3 Generator electrode	1 pc
11	Stirring bars($\varphi 7 \times 30$)	3 pcs
12	Drying cylinder	1 pc
13	Air tubing for waste(L=520mm, including joints)	1 pc
14	Air tubing for fill(L=580mm, including joints)	1 pc
15	Drain tube (L=800mm, including joints)	1 pc
16	Intake tube (L=800mm, including joints)	1 pc
17	Titration inlet tube (L=700mm, including joints)	1 pc
18	Pump tube (L=150mm, including joints)	1 pc
19	Cylinder for oven sampling tube	1 pc
20	Micro injector(100 μ L)	1 pc
21	Micro injector(1mL)	1 pc
22	Septum($\varphi 13 \times 3$)	20 pcs
23	Ejector tool (SC8.227.126)	1 pc
24	Titration outlet tube(L=600mm, including joints) and titration tip with anti-diffusion valve	1 pc
25	Burette device 10mL	1 pc

No.	Description	Quantity
26	Burette cover	1 pc
27	Fuse $\Phi 5 \times 20$ (4A)	3 pcs
28	Double magnetic ring power cord (black)	1 pc
29	M12-5 cable(Connect the Karl Fisher Titrator to the solvent manager)	1 pc
30	RS232-DB9 9 Pin male to male Serial Cross Cable-1 Meter	1 pc
31	USB communication cable-1.8 Meter	1 pc
32	USB to RS232 DB9 Serial Port Cable-1.5 Meter	1 pc
33	Touch pen (Capacitive screen)	1 pc
34	Instruction manual	1 pc
35	Sealing grease	1 pc
36	5A molecular sieve(500g)	1 pc