

Specifications

Model	DW-100MS	DW-100MS(plus)	DW-200MS	DW-200MS(plus)
Capacity	110g		210g	
Readability	0.001g			
Operating Temp	5°C-35°C			
Pan Size	Φ100mm			
Heating Source	Halogen Lamp			
Temp.Sensor	PT100			
Temp.Range	40°C-230°C			
Temp.Step	1°C			
Moisture Range	0.009%-100.0000			
Moisture Readability	0.01%			
Dry Range	100.00%-0.00%			
Dry Readability	0.01%			
Temp Setting	40°C-230°C By 1°C Step			
Time Setting	1-99min By Is			
Number of Storage	20			
Number Of History	20			
Heater	220V±15% 50HZ/110V±15% 60HZ/400W			
Balance	Input:220V 50HZ/110V 60HZ Output:9V			
Machine	Input:220V 50HZ/110V 60HZ Output:24V Dual-drive			
Host Dimension	215x195x415mm			
Indicator Dimension	200x75x135mm			
Operation Height	275mm			
The Chamber Heigh	35mm			
Display	7 inch touch panel			
Interface	RS232 (standard)		RJ45/USB (optional)	
Communication Mode	Cable	Wireless	Cable	Wireless
Li-Battery	No	10,000mAh	No	10,000mAh