

ICP 700T product shipping packaging configuration

Product Name : Inductively Coupled Plasma Optical Emission Spectrometer

Product Model : ICP700T

Number	Item name	Quantity	Unit	Remark	Confirm
Host					
1.1	ICP 700T Host	1	pc		
software					
2.1	ICP Spectrometer analysis software pieces	1	set	Lifetime free upgrades	
2.2	Product quality certificate	1	set	English	
2.3	Uplate	1	pc		
2.4	Warranty card	1	pc		
Accessories					
3.1	Network cable	1	pc	Category 5 twisted pair 3 rice	
3.2	plastic bucket	1	pcs	volume 5L	
3.3	Argon pressure reducing valve	2	pcs	Jierui RHM 152 IN -40	
3.4	PU hose	30	meter	10×6.5×1.75 mm	
3.5	PU hose	5	meter	6×4×1 mm	
3.6	PVC Transparent hose	0.5	meter	8 mm ×6 mm Transparent hose	
3.7	Plug power cord	1	pc	10A, 250V	
3.8	Grounding copper braided wire	15	meter	Cross-sectional area 4 mm ²	
3.9	brass tube	1	pc	long 500 mm , outer diameter 48 mm , wall thickness 3 mm	
3.10	filter cartridges	4	pcs	Outer diameter: 60 mm , inner diameter: 28 mm , length: 125 mm	
3.11	Torch	2	pcs	Quartz torch	
3.12	Double-tube spray chamber	2	pcs	Scott Spray chamber	
3.13	atomizer	1	pc	Imported atomizer	
3.14	High-frequency coil	1	pc		
3.15	Capillary injection tube	0.5	meter	D=2.38 mm , d=0.8 mm	
3.16	Peristaltic pump head hose	3	pcs	Φ1.65, long 45 cm , two cards	
3.17	Cartridge fuse	2	pcs	10A, 250V, 6×30 mm	
3.18	Cartridge fuse	2	pcs	6.3A, 250V, 5×20 mm	
3.19	Cylindrical spray chamber cover	2	pcs	polytetrafluoroethylene	
3.20	Silicone tube	2	meter	inner diameter 2 mm , outer diameter 4 mm , hardness 55	
3.21	Hose adapter	4	pcs	ZCJH -2*6 Straight through, inner diameter 6-2	
3.22	Hose adapter	4	pcs	ZC -1/16*3/32, inner diameter 1.6-2.4	
3.23	Ball joint clamp	1	pc	L00091181, Ø=10	

3.24	Quick Connectors	2	individual	Φ6	
3.25	Raw tape	1	roll	26 mm × 20 m, thickness: 0.1 to 0.2 mm	
3.26	O Type sealing ring	2	individual	32.5×3.55	
3.27	O Type sealing ring	4	individual	5.6×1.8	
3.28	Torch tooling	1	individual	self made	
3.29	cable ties	30	root	4×200	
Standard equipment					
4.1	Host	1	tower	Branded computers	
4.2	monitor	1	tower	Brand display	
4.3	printer	1	tower	Brand printer	
4.4	Cooling water tank	1	tower	Cooling capacity 2000W	
4.5	regulated power supply	1	tower	10kVA	
Packaging materials					
5.1	ICP 700T Packaging wooden boxes	1	set		
5.2	screw	5	individual	M5*8, hexagon socket, cross pan head	

