

ICP700T Pro Inductively Coupled Plasma Emission Spectrometer Configuration List

Serial Number	Accessory Name	Quantity	Unit	Specification And Model	Notes
I	Standard Configuration				
1	Host				
1.1	Inductively coupled plasma emission spectrometer	1	Set		
2	Software				
2.1	ICP spectrometer analysis software	1	Set	Free upgrade	
3	Accessory				
3.1	Network cable	1	Single	Category 5 twisted pair cables	
3.2	Plastic bucket	1	Piece	Volume 5L	
3.3	Argon pressure reducing valve	2	Piece	AR-8	
3.4	PU hose	30	Meter	Outside diameter 10mm	
3.5	PU hose	10	Meter	Outside diameter 6mm	
3.6	PVC transparent hose	1	Meter	Transparent hose	
3.7	Power cord	1	Single	10A	
3.8	USB extension cable	1	Single	Length 3m	
3.9	Grounding copper wire	15	Meter	4mm ²	
3.10	Ground rod	1	Piece	18*800	
3.11	Torch tube	2	Piece	1.2 quartz torch tube	
3.12	Atomizer	1	Piece	Quartz atomizer	
3.13	Double tube fog chamber	1	Piece	Scott fog chamber	
3.14	Cylindrical fog chamber cover	1	Piece	Polytetrafluoroethylene	
3.15	High-frequency coil component	1	Piece	Three turns silver plating	
3.16	Capillary sample tube	1	Meter	D=2.38mm, d=0.8mm	
3.17	Peristaltic pump head hose	3	Single	Length 45cm, Two clamps	
3.18	Tubular fuse	2	Piece	10A	
3.19	Tubular fuse	2	Piece	6.3A	
3.20	Silicone tube	4	Meter	Inside diameter 2mm, Outside diameter 4mm	
3.21	Hose adapter	4	Piece	Inside diameter 6-2	
3.22	Hose adapter	4	Piece	Inside diameter 1.6-2.4	
3.23	Spherical interface clamp	1	Piece	Stainless steel	

3.24	Quick coupling	3	Piece	Φ6	
3.25	Teflon tape	1	Roll	26mm×20m, Thickness: 0.1~0.2mm	
3.26	O-ring seal	4	Piece	32.5×3.55	
3.27	O-ring seal	8	Piece	5.6×1.8	
3.28	Torch pipe tooling	1	Piece	Cylinder	
3.29	Ribbon	50	Single	4×200	
3.30	Screw	5	Piece	M5, Cross dischead	
3.31	Screw	5	Piece	M4, Cross dischead	
3.32	T-shaped tee	2	Piece	10-10-10	
4	Standard Configuration				
4.1	Computer	1	Set	Brand	
4.2	Monitor	1	Set	Brand	
4.3	Printer	1	Set	Brand	
4.4	Cooling water tank	1	Set		
4.5	Regulated power supply	1	Set		
5	Delivery Attached Documenta tion				
5.1	U disk	1	Piece		
5.2	Product quality certificate	1	Part	English	
5.3	Warranty card	1	Part	English	