

ICP900T Inductively Coupled Plasma Emission Spectrometer Configuration List

Serial Number	Accessory Name	Quantity	Unit	Specification And Model	Notes
One	Standard Configuration				
1	Host				
1.1	Inductively coupled plasma emission spectrometer	1	Set		
2	Software				
2.1	ICP spectrometer analysis software	1	Set	Free upgrade	
3	Accessory				
3.1	Network cable	1	Piece	3m	
3.2	Plastic bucket	1	Piece	Volume 5L	
3.3	Argon pressure reducing valve	2	Piece	AR-8	
3.4	PU hose	30	Meter	10×6.5×1.75mm	
3.5	PU hose	5	Meter	6×4×1mm	
3.6	PVC transparent hose	1	Meter	8mm×6mm Transparent hose	
3.7	Power cord	1	Piece	10A, 250V	
3.8	USB extension cable	1	Piece	Length 3m	
3.9	Grounding copper wire	15	Meter	Sectional area 4mm ²	
3.10	Ground rod	1	Piece	18*800	
3.11	Torch tube	2	Piece	1.5 quartz torch tube	
3.12	Double tube fog chamber	1	Piece	Scott fog chamber	
3.13	Cylindrical fog chamber cover	1	Piece	Polytetrafluoroethylene	
3.14	Atomizer	1	Piece	Quartz atomizer	
3.15	Wavelength calibration solution	1	Bottle	1000ml	
3.16	High-frequency coil component	1	Piece	Three turns , silver plating	
3.17	Capillary sample tube	1	Meter	D=2.38mm, d=0.8mm	
3.18	Peristaltic pump head hose	3	Piece	Φ1.65, Length 45cm, Two clamps	
3.19	Tubular fuse	2	Piece	10A, 250V, 6×30mm	
3.20	Tubular fuse	2	Piece	6.3A, 250V, 5×20mm	
3.21	Silicone tube	4	Meter	Inside diameter 2mm, Outside diameter 4mm, H	

				ardness 55	
3.22	Hose adapter	4	Piece	ZCJH-2*6 direct connection, Inside diameter 6-2	
3.23	Hose adapter	4	Piece	ZC-1/16*3/32, Inside diameter 1.6-2.4	
3.24	Spherical interface clamp	1	Piece	L00091181, Ø=10	
3.25	Quick coupling	3	Piece	Φ6	
3.26	Teflon tape	1	Roll	Thickness: 0.1~0.2mm	
3.27	O-ring seal	4	Piece	32.5×3.55	
3.28	O-ring seal	8	Piece	5.6×1.8	
3.29	Torch pipe tooling	1	Piece	Cylinder	
3.30	Nylon cloth ventilation duct	10	Meter	Grey, Inside diameter 125mm	
3.31	Stainless steel pipe clamp	2	Piece	118-140mm	
3.32	Ribbon	50	Piece	4×200	
3.33	Screw	5	Piece	M5, Cross dischead	
3.34	Screw	5	Piece	M4, Cross dischead	
3.35	T-shaped tee	2	Piece	10-10-10	
3.36	Filter tank filter element	4	Piece	Outside diameter:60mm, Inside diameter:28mm, Length: 125mm	
4	Standard Configuration				
4.1	Computer	1	Set	Brand	
4.2	Monitor	1	Set	Brand	
4.3	Printer	1	Set	Brand	
4.4	Cooling water tank	1	Set		
4.5	Regulated power supply	1	Set		
5	Delivery Attached Documentation				
5.1	U disk	1	Piece		
5.2	Product quality certificate	1	Piece	English	
5.3	Warranty card	1	Piece	English	