

DWVP-30pro DWGM-4C

Corrosion Resistant Diaphragm Vacuum Pump



DWVP-30pro



DWGM-4C

Specifications

Model	DWVP-30pro	DWGM-4C
Maximum pumping speed at 50/60 Hz	2.1m³/h	3m³/h
Limit vacuum	5mbar / 3.6torr	
Number of pump heads/series	4 / 3	
Operating ambient temperature	10 - 40°C	
Storage ambient temperature	-10 - 60°C	
Maximum outlet absolute pressure	1.1bar	
Air inlet interface (DN)	11mm / 1/2 Hose joint	
Air outlet interface (DN)	11mm / 1/2 Hose joint	
Rated motor power	0.25kW	
Adjustable speed at 50/60 Hz	30 ~ 1800 RPM/min	30 ~ 2000 RPM/min
Speed of air extraction		1013 mbar: 50L/min 900 mbar: 45 L/min 750 mbar: 40 L/min 600 mbar: 35 L/min 480 mbar: 30 L/min 360 mbar: 25 L/min 250 mbar: 20 L/min 150 mbar: 15 L/min 50 mbar: 10 L/min 5 mbar: 5L/min
Degree of protection	IP 40	
Overall size	360*260*260mm	
Weight	16kg	